

PREFILLED SYRINGE KITS FOR VASCULAR ACCESS DEVICE FLUSH AND LOCK PROCEDURES

Each PraxiFlow™ kit consists of one Praxiject™ 0.9% NaCl syringe for Flush and one CitraFlow™ Sodium Citrate syringe for Lock of Vascular Access Devices. All the proven benefits of Praxiject™ and CitraFlow™ prefilled syringes in one convenient easy to use package.



- Easy to open package. Designed for use in a sterile field
- Ready to use prefilled syringes. Saves time and manipulations
- Each package contains one VAD Flush and one VAD Lock syringe
- Lower impact on environment. Less packaging than individually sold products.
- Available in 10ml syringe size for flushing as well as various Sodium Citrate concentrations for Locking

For further information, please contact your representative.

Catalog #	Description	Quantity/case
3215	1 Praxiject™ 0.9% NaCl syringe. (10ml in 10cc syringe) 1 CitraFlow™ 46.7% Sodium Citrate syringe. (3ml in 5cc syringe)	75 packages / cs
3225	1 Praxiject™ 0.9% NaCl syringe. (10ml in 10cc syringe) 1 CitraFlow™ 30% Sodium Citrate syringe. (3ml in 5cc syringe)	75 packages / cs
3255	1 Praxiject™ 0.9% NaCl syringe. (10ml in a 10cc syringe) 1 CitraFlow™ 4% Sodium Citrate syringe. (3ml in 5cc syringe)	75 packages / cs